

*All-in-one cleaning tool for KODAK MIN-R Screens and Cassettes*



# *Introducing* **KODAK MIN-R** Screen Cleaner Wipes



**K**odak Min-R screen cleaner wipes help make it easier than ever to keep Kodak Min-R screens and cassettes free of dust and dirt that can cause minus-density artifacts on film.

- Wipes are pre-moistened with Kodak Min-R screen cleaner—so there's no need for separate screen cleaner and wipes.
- Lint-free wipes significantly reduce concerns about artifacts caused by lint from cleaning cloths.
- A measured amount of screen cleaner in every wipe guards against screen staining or damage due to application of too much cleaner.



Another example of Kodak's continuing commitment to mammography customer satisfaction.

To include Kodak Min-R screen cleaner wipes in your next order, use CAT No. 528 0094. Or for more information, see your Kodak sales representative, call 1-800-328-2910, or visit [www.kodak/go/mammo](http://www.kodak/go/mammo).

**HEALTH IMAGING**  
A BETTER VIEW OF LIFE.



FIRST-CLASS MAIL  
US POSTAGE  
**PAID**  
ROCHESTER, NY  
PERMIT NO. 6

ADDRESS SERVICE REQUESTED