



**Kodak**

## Miniloader 2000/2000P

### Roomlight film handling in-house or on the road

The space-saving KODAK Miniloader 2000 and 2000P automatically load and unload KODAK MIN-R 2 Cassettes in full roomlight in just seconds. Now you can spend less time in the dark and more time with your patients.

#### Two configurations are available:

► **Standalone**

The KODAK Miniloader 2000, with its lightweight receive magazine and small footprint, is ideal for van installation when on-board film processing is not required.

► **Integrated**

A processor docked to the KODAK Miniloader 2000P allows automatic processing of exposed films for maximum convenience with easy image access and consistent image quality.

—If you already own a KODAK MIN-R Mammography Processor, upgrade with a KODAK Miniloader 2000P. It can also be equipped with a darkroom film feeder, and is suitable for van installation.

—The KODAK MIN-R Mammography Integrated Processor/Miniloader is a complete, newly integrated system that includes the MIN-R Mammography Processor and the Miniloader 2000P, as well as all the accessories necessary for a complete installation.

#### Lightweight receive magazine

Weighing only 1.4 kg (previously 4.2 kg), the lightweight receive magazine is designed to increase technologist comfort and satisfaction. In addition, the receive magazine improves workflow with a serial unload feature. High-quality images are assured with its light-tight design—no need to go to the darkroom to transfer films.

#### 2000P can serve a darkroom, too

Gain added cost savings and productivity by also processing films from an adjacent darkroom with optional through-the-wall feeding.

- Eliminate the need for a separate darkroom processor
- Simplify sensitometry procedures

#### User-friendly configuration

Front-accessible film-supply and receive magazines simplify workflow. An easy-to-read control panel displays operating and error conditions, as well as other key data. It even tracks film consumption on a daily, weekly, monthly, and annual basis.

#### The integrated approach to superior performance

Carestream Health's products, service, and support provide an integrated and comprehensive approach to meeting your needs with superior performance. From film processing to laser imaging, quality control to troubleshooting support, patient and technologist education to marketing expertise, Carestream Health leads the way. In the fight against breast cancer and as an advocate for women's health, you won't find a better partner.

# Kodak Miniloader 2000/2000P



## Specifications

### Cycle time

12 to 15 sec. from cassette insertion to reload, depending on film size and magazine position

### Dimensions

#### 2000 configuration for mobile units

Width: 55 cm +/- 1 cm /21.6 in.

Depth: 60 cm +/- 1 cm /23.6 in.

Height: 124 cm/49 in.

Weight: 220 kg/484 lb

#### 2000P configuration with integrated processor

Width: 73 cm/28.3 in.

Depth: 78 cm/31.5 in.

Height: 138 cm/55 in.

Weight: 250 kg/550 lb (w/o processor)

### Electrical

110/240 V, single phase

300 W max. consumption

UL listed; CSA, TUV-GS, and CE certified

### Cassettes

Accepts KODAK MIN-R 2 cassettes (18 x 24 cm and 24 x 30 cm).

## Ordering information

### CONFIGURATION

### CAT NO.

#### Standalone

- ▶ KODAK Miniloader 2000 . . . . 721 8969

#### Integrated

- ▶ KODAK Miniloader 2000P . . . 721 8977
- ▶ KODAK MIN-R Integrated Processor/Miniloader . . . . . 132 3526

## 2000 accessories

- ▶ KODAK Miniloader 2000 humidifier, model 2
- ▶ KODAK Miniloader 2000 van kit, model 2
- ▶ KODAK Miniloader 2000 video film holder kit
- ▶ KODAK Miniloader 2000 lightweight receive magazine
- ▶ KODAK Miniloader 2000 supply magazine

## 2000P accessories

- ▶ KODAK Miniloader 2000P serial unload adapter kit/for multiloaders 300/700/7000
- ▶ KODAK Miniloader 2000P van kit
- ▶ KODAK Miniloader 2000P darkroom film feeder
- ▶ Mounting stand for KODAK MIN-R Mammography Processor
- ▶ KODAK Miniloader 2000P/MIN-R Processor\* ready kit
- ▶ KODAK Miniloader 2000P humidifier
- ▶ KODAK Miniloader 2000P supply magazine

Note: Processor and accessories are sold separately.

\*KODAK MIN-R Mammography processor replaces the M35-A/M processor.

Feature	KODAK Miniloader 2000	KODAK Miniloader 2000P
Lightweight receive magazine	Yes	No
Film processor docking capabilities	No	Yes
Load-only function for loading empty cassette	Yes	Yes
Unload-only function to empty cassette for cleaning or inspection	Yes	Yes
Serial unloading for rapid unloading of sequential films from receive magazine	No	Yes
Accepts 8 x 10-in. KODAK X-OMATIC video film holders for ultrasound applications	No	Yes

## More information

For details on the KODAK Miniloader 2000/2000P, or to order, contact your Carestream Health representative, or call 1-877-865-6325, ext. 626.

➔ [www.carestreamhealth.com](http://www.carestreamhealth.com)

**Kodak**

Licensed Product

Carestream Health, Inc.  
150 Verona Street  
Rochester, NY 14608

Carestream Health Canada Company  
6 Monogram Place  
Suite 200  
Toronto, Ontario M9R 0A1  
CANADA

**Carestream**   
HEALTH